



# Recombinant Lactococcus lactis subsp. cremoris Enolase 1 (eno1)

<b>Product Code</b>	CSB-EP312687LBE
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q032H8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Lactococcus lactis subsp. cremoris (strain SK11)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTVTIENIHA REIFDSRGNP TVEVDVRLTD GTLGRAAVPS GASTGDREAV ELRDGGARLQ GKDVSKAVAN VNGEIYEALK GQSPFNQAKL DHLMIELDGT KNKSRLGANA ILGVSMAIAR AAANSEKIPL YRYLGGVDLE LPQPPFNVIN GGVHADSGID VQEFLITPVK RASFRDGLEK IANIHHTLKK ILADKGLETA VGDEGGFAPK LGSTENAIAT LYQAIERAGY VPGEEIAIAI DPASSEFYDD KEKVYHFEGQ KL TSAELLY YEGLVEKYPA LISIEDGFSE HDWAGFAAQT KVQGGKIQLV GDDIFVTNPE IFKEGIQKGV ANAILIKLNQ IGTVTEAIEA ISLARKAGYK TMISHRSGET VDSFIADFAV AMHAGQIKTG SMARSERVEK YNQFLRIEEE LLGLEVTK
<b>Source</b>	E.coli
<b>Target Names</b>	eno1
<b>Protein Names</b>	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
<b>Expression Region</b>	1-418
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.