



Recombinant *Lactococcus lactis* subsp. *cremoris* Enolase 1 (eno1)

Product Code	CSB-YP312687LBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q032H8
Product Type	Recombinant Protein
Immunogen Species	<i>Lactococcus lactis</i> subsp. <i>cremoris</i> (strain SK11)
Purity	>85% (SDS-PAGE)
Sequence	MTVTIENIHA REIFDSRGNP TVEVDVRLTD GTLGRAAVPS GASTGDREAV ELRDGGARLQ GKDVSKAVAN VNGEIYEALK GQSPFNQAKL DHLMIELDGT KNKSRLGANA ILGVSMAIAR AAANSEKIPL YRYLGGVDLE LPQPFNVIN GGVHADSGID VQEFLITPVK RASFRDGLEK IANIHHTLKK ILADKGLETA VGDEGGFAPK LGSTENAIAT LYQAIERAGY VPGEEIAIAI DPASSEFYDD KEKVYHFEGQ KL TSAELLY YEGLVEKYPA LISIEDGFSE HDWAGFAAQT KVQGGKIQLV GDDIFVTNPE IFKEGIQKGV ANAILIKLNQ IGTVTEAIEA ISLARKAGYK TMISHRSGET VDSFIADFAV AMHAGQIKTG SMARSERVEK YNQFLRIEEE LLGLEVTK
Source	Yeast
Target Names	eno1
Protein Names	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.