



Recombinant Mouse Proto-oncogene Wnt-1 (Wnt1)

Product Code	CSB-EP026128MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04426
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>ANS SGRWWGIVNI ASSTNLLTDS KSLQLVLEPS LQLLSRKQRR LIRQNPGLIH SVSGGLQSAV RECKWQFRNR RWNCPTAPGP HLFKIVNRG CRETAFIFAI TSAGVTHSVA RSCSEGSIES CTCDYRRRGP GPDWHWGGC SDNIDFGRLF GREFVDSGEEK GRDLRFLMNL HNNEAGRRTV FSEMRQECKC HGMSGSCTVR TCWMRLPTLR AVGDVLRDRF DGASRVLYGN RGSNRASRAE LLRLEPEDPA HKPPSPHDLV YFEKSPNFCT YSGRLGTAGT AGRACNSSSP ALDGCELLCC GRGHRTRTQR VTERCNC TFH WCCHVSCRNC THTRVLHECL</p>
Source	E.coli
Target Names	Wnt1
Protein Names	Recommended name: Proto-oncogene Wnt-1 Alternative name(s): Proto-oncogene Int-1
Expression Region	28-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	<p>The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It is very conserved in evolution, and This protein is known to be 98% identical to the mouse Wnt1 protein at the amino acid level. The studies in mouse indicate that the Wnt1 protein functions in the induction of the mesencephalon and cerebellum. This gene was originally considered as a candidate gene for Joubert syndrome, an autosomal recessive disorder with cerebellar hypoplasia as a leading feature. However, further studies suggested that the gene mutations might not have a significant role in Joubert syndrome. This gene is clustered with another family member, WNT10B, in the chromosome 12q13 region.</p>
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.