



# Recombinant Rat Interleukin-3 (IL3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP011652RA  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P04823  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | SDR GSDAHLLRT LDCRTIALEI LVKLPVSGLN NSDDKANLRN<br>STLRRVNLDE FLKSQEEFDS QDTTDIKSKL QKLKCCIPAA ASDSVLPGVY<br>NKDLDDFKKK LRFYVIHLKD LQPVSVSRPP QPTSSSDNFR PMTVEC  |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | IL3   |
| <b>Protein Names</b>     | Recommended name: Interleukin-3 Short name= IL-3 Alternative name(s):<br>Hematopoietic growth factor Mast cell growth factor Short name= MCGF<br>Multipotential colony-stimulating factor P-cell-stimulating factor   |
| <b>Expression Region</b> | 28-166  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | This protein is a potent growth promoting cytokine. This cytokine is capable of supporting the proliferation of a broad range of hematopoietic cell types. It is involved in a variety of cell activities such as cell growth, differentiation and apoptosis. This cytokine has been shown to also possess neurotrophic activity, and it may be associated with neurologic disorders.                       |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |