



# Recombinant Rat Proliferating cell nuclear antigen (Pcna)

<b>Product Code</b>	CSB-BP017621RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04961
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFEARLIQGS ILKKVLEALK DLINEACWDI SSGGVNLQSM DSSHVSLVQL TLRSEGFDTY RCDRNLAMGV NLTSMKILK CAGNEDIITL RAEDNADTLA LVFEAPNQEK VSDYEMKLM DLDVEQLGIPE QEYSCVVKMP SGEFARICRD LSHIGDAVVI SCAKDGVKFS ASGELGNGNI KLSQTSNVDK EEEAVSIEMN EPVQLTFALR YLNFFTKATP LSPTVTLSMS ADVPLVVEYK IADMGHLKYY LAPKIEDEEG S
<b>Source</b>	Baculovirus
<b>Target Names</b>	Pcna
<b>Protein Names</b>	Recommended name: Proliferating cell nuclear antigen Short name= PCNA Alternative name(s): Cyclin
<b>Expression Region</b>	1-261
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein is found in the nucleus and is a cofactor of DNA polymerase delta. The encoded protein acts as a homotrimer and helps increase the processivity of leading strand synthesis during DNA replication. In response to DNA damage, this protein is ubiquitinated and is involved in the RAD6-dependent DNA repair pathway. Two transcript variants encoding the same protein have been found for this gene. Pseudogenes of this gene have been described on chromosome 4 and on the X chromosome.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.