



Recombinant Human Interferon alpha-6 (IFNA6)

Product Code	CSB-MP011042HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P05013
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SLDCDLPQTH SLGHRRTMML LAQMRRISLF SCLKDRHDFR FPQEEDGNQ FQKAEAISVL HEVIQQTFLN FSTKDSSVAW DERLLDKLYT ELYQQQLNDLE ACVMQEVWVG GTPLMNEDSI LAVRKYFQRI TLYLTEKKYS PCAWEVVRAE IMRSFSSSRN LQERLRRKE
Source	Mammalian cell
Target Names	IFNA6
Protein Names	Recommended name: Interferon alpha-6 Short name= IFN-alpha-6 Alternative name(s): Interferon alpha-54 Interferon alpha-K Short name= LeIF K
Expression Region	21-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.