



# Recombinant Pseudomonas putida Phthalate 4,5-dioxygenase oxygenase subunit (pht3)

<b>Product Code</b>	CSB-BP313137FFZ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q05183
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pseudomonas putida (Arthrobacter siderocapsulatus)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLTPEENLLL CRVEGDAPMG QMMRRHWTPV CLLEEVSEPD GTPVRARLFG EDLVVFRDTD GRVGVMD EYC PHRRVSLIYG RNENGLRCL YHGWMKMDVDG NVVEMVSEPA ASNMCQKVKH TAYKTREWGG FVWAYMGPQD AIPEFVPPAW APHEHVRVSI AKAIIPC NWA QILEGAIDSA HSSSLHSSDF VPARVGGAEA TSKNWL RPST DKAPRMQVER TSYGFRYAAL RRPIQNAATS EYVRSTVFVA PATALIPPNN LYNVANINVP IDDTHTAFYF MAWGNPDNTP ETETWRKFLG QQVGIDLDDS YRPLRNDG NR FFQDREAMKN GNFTGIKGFP NQDIAMWV TM GPIADRSDER LGASDLAVVE FRRVMLDALA AFQAGESAIG TGEKAIPSRI CSFQAIVSKD IDWRDYQARY VWALDDANIV AEPDYE VHT
<b>Source</b>	Baculovirus
<b>Target Names</b>	pht3
<b>Protein Names</b>	Recommended name: Phthalate 4,5-dioxygenase oxygenase subunit EC=1.14.12.7
<b>Expression Region</b>	1-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.