



# Recombinant *Saccharomyces cerevisiae* Protein PAL1 (PAL1)

<b>Product Code</b>	CSB-BP313183SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q05518
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MENRNSSTSS RPFVNNPFR NATVDSSINQ YKNDSQFQEW AKNQSRNSF DMPQLNTRTS SQLSFPNPE DEPQRNADQQ GAFYSGLESF SSGSLPPSR PLSSKNPFLD DVSSATDFRR SPPPVS RNKN HPTAKEEKEQ LRQRYLEESD VSTVGNTREN TDLPPSYEEI TSTNGSRRAY PKEKVS RPSS HREHSNSGTY ISRRSSSHHH REASSSSTPS KKGKRKSKVI VPKNVD TIDK LDVTGLFGGS FHHDGPFDAV TPHRNKNNKA APVLAFPVDG PNSTIGGAST KKSALDEVFG RDDTDDSDIY QYSSQTLRRG GDTQDAIKAN VGNVQQMDAK NKTEL VHGPV TAGLGSSTFL DGAPASSAAI RNDIKAHSYH NRNGGLQRNK SLSQRLGLGG SGDSNAPMTG VRRNLSLSRD NYDVGH SNEG VRRSKTVNSP NRTHKS NYTT DFDGQDDHNE DEEDVYLGVR YNEPNMKKKS TGSKLLSRVK SLKVGRKSQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	PAL1
<b>Protein Names</b>	Recommended name: Protein PAL1 Alternative name(s): Pears and lemons protein 1
<b>Expression Region</b>	1-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.