



Recombinant *Saccharomyces cerevisiae* Protein PAL1 (PAL1)

Product Code	CSB-YP313183SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05518
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MENRNSSTSS RPFVNNPFR NATVDSSINQ YKNDSQFQEW AKNQSRNSF DMPQLNTRTS SQLSFPNIPE DEPQRNADQQ GAFYSGLESF SSGSLPPSR PLSSKNPFLD DVSSATDFRR SPPPVS RNKN HPTAKEEKEQ LRQRYLEESD VSTVGNTREN TDLPPSYEEI TSTNGSRRAY PKEKVS RPSS HREHSNSGTY ISRRSSSHHH REASSSSTPS KKGKRKSKVI VPKNVDTIDK LDVTGLFGGS FHHDGPFDAV TPHRNKNNKA APVLAFPVDG PNSTIGGAST KKSALDEVFG RDDTDDSDIY QYSSQTLRRG GDTQDAIKAN VGNVQQMDAK NKTEL VHGPV TAGLGSSTFL DGAPASSAAI RNDIKAHSYH NRNGGLQRNK SLSQRLGLGG SGDSNAPMTG VRRNLSLSRD NYDVGH SNEG VRRSKTVNSP NRTHKSNYTT DFDGQDDHNE DEEDVYLGVR YNEPNMKKKS TGSKLLSRVK SLKVGRKSQ
Source	Yeast
Target Names	PAL1
Protein Names	Recommended name: Protein PAL1 Alternative name(s): Pears and lemons protein 1
Expression Region	1-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.