



Recombinant Salmonella typhimurium Putative cobalt-precorrin-6A synthase [deacetylating] (cbiD)

Product Code	CSB-MP313199SXB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05628
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	MSELSFDAPV WHHGKALRKG YTTGSCATAA AKVAALMVLR QHLIHQVSIV TPSGVTLCLN VESPHIEGQQ AIAAIRKDGG DDVDATHGML IFARVTLNDS GEITLTGGEG IGTVTRKGIG LPLGSAAINR TPRHTIESAV REAIGPARGA DVEIFAPEGE ARAQKTYNSR LGILGGISII GTTGIVTPMS EESWKRSLSL ELEIKRASGL TRVILVPGNH GERFVREQMG VDTQAVVTMS NFVGYMIEEA VRLGFCQIVL VGHPGKLIKI AAGIFHTHSH IADARMETLV AHLALLGAPL ELLTLVSDCD TTEAAMEHIE AYGFGHIYNH LARRICLRVM QMLRFTKTPP VCDAILFSFD NHILGSNRPV DEIAKELQC
Source	Mammalian cell
Target Names	cbiD
Protein Names	Recommended name: Putative cobalt-precorrin-6A synthase [deacetylating] EC= 2.1.1.-
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.