



# Recombinant *Saccharomyces cerevisiae* RNA-binding suppressor of PAS kinase protein 1 (RBS1)

<b>Product Code</b>	CSB-YP313202SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q05672
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTALDSRNWG LTPAMETGLF QKPQDRIFII ELENSIVSFI NSNTESFQLR PMNSYYRLLS HQIAEYHNLN HVLARTQDSC VILFKGENFQ KIEGKPLLQE LQLNKKPEEC ASSSESIEKS NNNRIFRILK RKEVGNEDY KIDGNTRTPN SNLTANSNKD QKIEIDDKSS TDLEQERIEK ERLYEQRKQE IFDKLNKSED DVKSTNSSGS SDSDNEWSDW LNGDDSNTQT SNGSVSSSSP FNSSVTTIQV NKPQQQFYDS RRGRGGRRRG TNNYKDAYRG QSRRNKENG YQSGYSSPYL VYPPPQMGGN SLPTYPLMYN PAGPAPGPAP SPMVMGNNTV FMNPYMYNMN PQGSCSFGTP IPMYPPYQYQ YQYQYNTQYH SGPYSNTPSY NSNNYTRSSA NKYHHFQGKN SYSGAIPKRS DDSNSNKNEG IRRASVEGSP SSRDTSVEM KFDKLN
<b>Source</b>	Yeast
<b>Target Names</b>	RBS1
<b>Protein Names</b>	Recommended name: RNA-binding suppressor of PAS kinase protein 1
<b>Expression Region</b>	1-457
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.