



Recombinant Outer membrane protein IIIA (ropA)

Product Code	CSB-EP313222RKR-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05811
Product Type	Recombinant Protein
Immunogen Species	Rhizobium leguminosarum bv. viciae
Purity	>85% (SDS-PAGE)
Sequence	ADAIVAAE PEPVEYVRVC DAYGTGYFYI PGTETCLKIE GYIRFQVNVG DNPGGDNDSD WDAVTAVRFS SRKSDTEYGP LTGVIVMQFN ADNASDQDAI LDSAYLDVAG FRAGLFYSWW DDGLSGETDD IGSVVTLHNS IRYQYESGTF YAGLSVDELE DGVYQGTFTP GVIPGTTDFT ADDGPNNVGV AFGIGGTAGA FSYQVTGGWD VDNEGAIIRA MGTVEIGPGT FGLAGVYSSG PNSYYSSAEW AVAAEYAIKA TDKLKITPGR WHGHVPEDFD GLGDAWKVGL TVDYQIVENF YAKASVQYLD PQDGEDSTSG YFACSVRSNH LVDAPGLRIG STTISF
Source	E.coli
Target Names	ropA
Protein Names	Recommended name: Outer membrane protein IIIA Short name= OMP IIIA
Expression Region	23-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.