



Recombinant *Saccharomyces cerevisiae* Central kinetochore subunit NKP2 (NKP2)

Product Code	CSB-YP313269SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06162
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNSEQLLHNY VSDSLTTLI SFQEFKQQLQ SYTSDEQQLQ HWYELLQARD ARVTSELEAR IKQFFITLRS RLLRFLESEQ LSHLSLETI IDALYKINDL LQQLRLQILDD AIQEKTSELA EFENMVRSPS AGDNAIPGLL QIIQSYINLL EEN
Source	Yeast
Target Names	NKP2
Protein Names	Recommended name: Central kinetochore subunit NKP2 Alternative name(s): Non-essential kinetochore protein 2
Expression Region	1-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.