



Recombinant *Saccharomyces cerevisiae* Putative hexokinase YLR446W (YLR446W)

Product Code	CSB-MP313275SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06204
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTIESTLARE LESLILPADS IVNVVDQFQE ELLSRLQTNT ISMLPQCLVP DKRSRWNPED KILTIDFGGT RLKFAIISLP QIVIEYNDAF ELTYNIVDSN FFNQIYTIC TRLAANGYIK KKNESSEASK FFVSVTFSFP LNPEGEVVAM GKGFVMTDTL QGSTVKQLIQ SSFHRIISEN IEEFFCTMNV CHVINDAIAV SLTSKFICEN DSISLIIGTG TNACFEVPG YLPPFKRDAL RETLPSSYNK ETLNFKHVLI NSEIGFIGKN VIALQPFDIH GAISYEMPLE CVTSGKWLPL SLKNILLQYN IIPKNFPVEF NGELVCQLAE DCTNAWFENE HYALICQIAR LLIKRAAFYV AAIVQAIDII TGCKNYNFIH IGYVGSFLHN SNFYREQIKY YSSIIHKLQF LNHSNLLGAA IATYLNKSDN QVQ
Source	Mammalian cell
Target Names	YLR446W
Protein Names	Recommended name: Putative hexokinase YLR446W EC= 2.7.1.1
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.