



Recombinant *Saccharomyces cerevisiae* F-box protein YDR306C (YDR306C)

Product Code	CSB-EP313308SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06640
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MANKSRPKKI KAPYRKYVAG EGFSSSTRNDN KAKEFTITIP EDAELIETPQ GSYYYDETND TIVKLTRLSN EKKDKKGRKQ SPSSSSTSSS KGEKNGKVIE SEEARMHSVS VKMVLPWEIQ HRIIHYLDIP EKEEKLNKTA NGKKTGGIN MNYLLVCRNW YAMCLPKLYY APALTSKNFN GFVDTIINK KKNLGHYVFE LNLSTILQSG RNSFVSKLLR RCCSNLTKFI APQTSFGYAP LISLKSCHDL KFLDLGLVSE TVKLKELFSA IKNFTKLTHL SFPRSSIDCQ GFQDIQWPQN LRYLKLSSGI TNEFVIDTKW PTTITTLEFS YCPQITELSI YSLLSQIGDN LKHLFFHYPM PSLAENSLDH VFTYCANLIS LQLMVDYCSK WCFSEFMLSK LVEYDRPLKT LYLECSGSLG LASKIHPDDL TIAILESRLP CLKNICVSPK LGWNMKSDEV ADLVVSLEDQ DGSLYLNY</p>
Source	E.coli
Target Names	YDR306C
Protein Names	Recommended name: F-box protein YDR306C
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.