



Recombinant Succinyl-CoA ligase [ADP-forming] subunit alpha (sucD)

Product Code	CSB-MP313333EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGF0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	SILIDKNTK VICQGFTGSQ GTFHSEQAIA YGTMVGGVT PGKGGTTHLG LPVFNTVREA VAATGATASV IYVPAPFCKD SILEAIDAGI KLIITITEGI PTLDMLTVKV KLDEAGVRMI GPNCPGVITP GECKIGIQPG HHHKPGKVGI VSRSGTLTYE AVKQTTDYGF GQSTCVGIGG DPIPGSNFID ILEMFEKDPQ TEAIVMIGEI GGSAAAAAAA YIKEHVTKPV VGYIAGVTAP KGKRMGHAGA IIAGGKGTAD EKFAALEAAG VKTVRSLADI GEALKTVLK
Source	Mammalian cell
Target Names	sucD
Protein Names	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit alpha EC=6.2.1.5 Alternative name(s): Succinyl-CoA synthetase subunit alpha Short name= SCS-alpha
Expression Region	2-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.