



Recombinant Protein sirB1 (sirB1)

Product Code	CSB-MP313344EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGM6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MRSLADFEFN KAPLCEGMIL ACEAIRRDFP SQDVYDELER LVSLAKEEIS QLLPLEEQLE KLIALFYGDW GFKASRGVYR LSDALWLDQV LKNRQGS AVS LGAVLLWVAN RLDLPLLPVI FPTQLILRIE CPDGEIWLIN PFNGESLSEH MLDVWLKGN I SPSAELFYED LDEADNIEVI RKLLDTLKAS LMEENQMELA LRTSEALLQF NPEDPYEIRD RGLIYAQLDC EHVALNDLSY FVEQC PEDPI SEMIRAQINN IAHKHIVLH
Source	Mammalian cell
Target Names	ychA
Protein Names	Recommended name: Protein sirB1
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.