



Recombinant Protein W (P/V/C)

Product Code	CSB-EP313495NDT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C1C7
Product Type	Recombinant Protein
Immunogen Species	Nipah virus
Purity	>85% (SDS-PAGE)
Sequence	MDKLELVNDG LNIIDFIQKN QKEIQKTYGR SSIQQPSIKD QTKAWEDFLQ CTSGESEQVE GGMSKDDGDV ERRNLEDLSS TSPTDGTIGK RVSNTRDWAE GSDDIQLDPV VTDVVYHDHG GECTGYGFTS SPERGWSDYT SGANNGNVCL VSDAKMLSYA PEIAVSKEDR ETDLVHLENK LSTTGLNPTA VPFTLRNLSD PAKDSPVIAE HYYGLGVKEQ NVGPQTSRNV NLDSIKLYTS DDEEADQLEF EDEFAGSSSE VIVGISPEDE EPSSVGGKPN ESIGRTIEGQ SIRDNLQAKD NKSTDVPGAG PKDSAVKEEP PQKRLPMLAE EFECGSEDPE IIRELLKENS LINCQQGKDA QPPYHWSIER SISPKTEIV NGAVQTADRQ RPGTPMPKSR GIPIKKGAT RNIHLLGRKT CLGRRVVQPG MFEDHPPTKK ARVSMRMSN
Source	E.coli
Target Names	P/V/C
Protein Names	Recommended name: Protein W
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.