



Recombinant Human Tetranectin (CLEC3B)

Product Code	CSB-BP005531HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P05452
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	EPPTQKPKK IVNAKKDVVN TKMFEELKSR LDTLAQEVAL LKEQQALQTV CLKGTKVHMK CFLAFTQTKT FHEASEDCIS RGGTLGTPQT GSENDALYEY LRQSVGNEAE IWLGLNDMAA EGTWVDMTGA RIAYKNWETE ITAQPDDGGKT ENCAVLSGAA NGKWFDKRCR DQLPYICQFG IV
Source	Baculovirus
Target Names	CLEC3B
Protein Names	Recommended name: Tetranectin Short name= TN Alternative name(s): C-type lectin domain family 3 member B Plasminogen kringle 4-binding protein
Expression Region	22-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.