



Recombinant *Oryza sativa* subsp. japonica Serine/threonine-protein kinase SAPK3 (SAPK3)

Product Code	CSB-EP313613OFG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C5D6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEERYEALKE LGAGNFGVAR LVRDKRSKEL VAVKYIERGK KIDENVQREI INHRS LRHPN IIRFKEVCLT PTHLAIVMEY AAGGELFEQI CTAGRFSEDE ARYFFQQLIS GVSYSCHSLEI CHRDLKLENT LLDGSPTPRV KICDFGYSKS ALLHSKPKST VGTPAYIAPE VLSRKEYDYGK VADVWSCGVT LYVMLVGSYP FEDPGDPRNF RKTISRILGV QYSIPDYVRV SSDCRLLSQ IFVADPSKRI TIPEIKKHTW FLKNLPKEIS EREKADYKDT DAAPPTQAVE EIMRIIQEAK VPGDMAAADP ALLAELAELK SDDEEEAADE YDTY
Source	E.coli
Target Names	SAPK3
Protein Names	Recommended name: Serine/threonine-protein kinase SAPK3 EC= 2.7.11.1 Alternative name(s): Osmotic stress/abscisic acid-activated protein kinase 3 Protein kinase REK
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.