



Recombinant *Oryza sativa* subsp. *indica* Prolamin PPROL 14E

Product Code	CSB-EP313615OFF-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C5E5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	Q FDVLGQSYRQ YQLQSPVLLQ QQVLSPYNEF VRQQYGIAAS PFLQSAAFQL RNNQVWQQLA LVAQQSHYQD INIVQAIQQ LQLQQFGDLY FDRNLAQAQA LLAFNVPSRY GIYPRYYGAP STITTLGGVL
Source	E.coli
Protein Names	Recommended name: Prolamin PPROL 14E Alternative name(s): Prolamin PPROL 14 Prolamin PPROL 4A Prolamin PPROL 7
Expression Region	20-150
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.