




# Recombinant Hepatitis delta virus genotype I Large delta antigen

<b>Product Code</b>	CSB-YP313702HFN
<b>Relevance</b>	Following virus entry into host cell, provides nuclear import of HDV RNPs thanks to its nuclear localization signal. Needs co-infection with hepatitis B virus to provide surface proteins, otherwise there is no packaging or budding. Packages the HDV ribonucleoprotein in hepatitis B virus empty particles. Interacts with both HDV genomic RNA and cytoplasmic tail of HBsAg. May inhibit viral RNA replication
<b>Abbreviation</b>	Recombinant Hepatitis delta virus genotype I Large delta antigen protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0C6L6
<b>Alias</b>	p27
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Hepatitis delta virus genotype I (isolate Italian) (HDV)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MSRSESRKNRGGREEILEQWVAGRKKLEELERDLRKTKKKLLKIEDENPWLGNIKGILGKKDKDGEGAPPAKRARTDQMEVDSGPRKRPLRGGFTDKERQDHRRKALENKKKQLSAGGKNLSKEEEEEELRRLTEEDERRRRVAGPPVGGVIPLEGGSRGAPGGGFVPSLQGVPEPFSRTGEGLDIRGNRGFPWDILFPADPPFSPQSC
<b>Research Area</b>	others
<b>Source</b>	Yeast
<b>Protein Names</b>	Recommended name: Large delta antigen Short name= L-HDAg Alternative name(s): p27
<b>Expression Region</b>	1-211aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	27.7kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.