



Recombinant Ureaplasma parvum serovar 3 Urease subunit beta (ureB)

Product Code	CSB-EP313760UAO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C7K8
Product Type	Recombinant Protein
Immunogen Species	Ureaplasma parvum serovar 3 (strain ATCC 700970)
Purity	>85% (SDS-PAGE)
Sequence	SGSSSQFSP GKLVPAINF ASGEIVMNEG REAKVISIKN TGDRPIQVGS HFHLFEVNSA LVFFDEKGNE DKERKVAYGR RFDIPSGTAI RFEPGDKKEV SIIDLAGTRE VWGVNGLVNG KLKK
Source	E.coli
Target Names	ureB
Protein Names	Recommended name: Urease subunit beta EC= 3.5.1.5 Alternative name(s): Urea amidohydrolase subunit beta
Expression Region	2-124
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.