



Recombinant Signal recognition particle protein (ffh)

Product Code	CSB-EP313780EGX-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGD8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MFDNLTDRLS RTLRNISGRG RLTEDNVKDT LREVRMALLE ADVALPVVRE FINRVKEKAV GHEVNKSLTP GQEFVKIVRN ELVAAMGEEN QTLNLAAQPP AVVLMAGLQG AGKTTSVGKL GKFLREKHKK KVLVVSADV RPAAIKQLET LAEQVGVDFF PSDVGQKQVD IVNAALKEAK LKFYDVLLVD TAGRLHVDEA MMDEIKQVHA SINPVETLFV VDMGTQDAA NTAKAFNEAL PLTGVVLTKV DGDARGGAAL SIRHITGKPI KFLGVGEKTE ALEPFHPDRI ASRILGMGDV LSLIEDIESK VDRAQAEKLA SKLKKGDGFD LNDFLEQLRQ MKNMGGMASL MGKLPGMGQI PDNVKSQMDD KVLVRMEAI NSMTMKERAK PEIIKGSRKR RIAAGCGMQV QDVNRLKQF DDMQRMMKKM KKGGMAMMMR SMKGMMPPGF PGR
Source	E.coli
Target Names	ffh
Protein Names	Recommended name: Signal recognition particle protein Alternative name(s): Fifty-four homolog
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.