



# Recombinant Streptococcus pyogenes Arginine deiminase (arcA)

<b>Product Code</b>	CSB-BP313805SMR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0C0B3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Streptococcus pyogenes
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TAQTPIHVY SEIGKLLKVL LHRPGKEIEN LMPDYLERLL FDDIPFLEDA QKEHDAFAQA LRDEGIEVLY LETLAAESLV TPEIREAFID EYLSEANIRG RATKKAIREL LMAIEDNQEL IEKTMAGVQK SELPEIPASE KGLTDLVESN YPFAIDPMPN LYFTRDPFAT IGTGVSLNHM FSETRNRETL YGKYIFTHHP IYGGGKVPMV YDRNETTRIE GGDELVLSKD VLAVGISQRT DAASIEKLLV NIFKQNLGFK KVLAFEFANN RKFMHLDTVF TMVDYDKFTI HPEIEGDLRV YSVTYDNEEL HIVEEKGLDA ELLAANLGVE KVDLIRCGGD NLVAAGREQW NDGSNTLTIA PGVVVVYNRN TITNAILESK GLKLIKIHGS ELVRGRGGPR CMSMPFERED I
<b>Source</b>	Baculovirus
<b>Target Names</b>	arcA
<b>Protein Names</b>	Recommended name: Arginine deiminase Short name= ADI EC= 3.5.3.6 Alternative name(s): Arginine dihydrolase Short name= AD Streptococcal acid glycoprotein
<b>Expression Region</b>	2-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.