



Recombinant Streptococcus pyogenes Arginine deiminase (arcA)

Product Code	CSB-YP313805SMR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0B3
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes
Purity	>85% (SDS-PAGE)
Sequence	TAQTPIHVY SEIGKLLKVL LHRPGKEIEN LMPDYLERLL FDDIPFLEDA QKEHDAFAQA LRDEGIEVLY LETLAAESLV TPEIREAFID EYLSEANIRG RATKKAIREL LMAIEDNQEL IEKTMAGVQK SELPEIPASE KGLTDLVESN YFPAIDPMPN LYFTRDPFAT IGTGVSLNHM FSETRNRETL YGKYIFTHHP IYGGGKVPMV YDRNETTRIE GGDELVLSKD VLAVGISQRT DAASIEKLLV NIFKQNLGFK KVLAFEFANN RKFMHLDTVF TMVDYDKFTI HPEIEGDLRV YSVTYDNEEL HIVEEKGLA ELLAANLGVE KVDLIRCGGD NLVAAGREQW NDGSNTLTIA PGVVVVYNRN TITNAILESK GLKLIKHGS ELVRGRGGPR CMSMPFERED I
Source	Yeast
Target Names	arcA
Protein Names	Recommended name: Arginine deiminase Short name= ADI EC= 3.5.3.6Alternative name(s): Arginine dihydrolase Short name= AD Streptococcal acid glycoprotein
Expression Region	2-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.