



Recombinant Protein rdoA (rdoA)

Product Code	CSB-BP313812EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C0K4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNNSAFTFQT LHPDTIMDAL FEHGIRVDSG LTPLNSYENR VYQFQDEDRR RFVVKFYRPE RWTADQILEE HQFALQLVND EVPVAAPVAF NGQTL LN HQG FYFAVFP SVG GRQFEADNID QMEAVGRYLG RMHQTGRKQL FIHRPTIGLN EYLIEPRKLF EDATLIPSGL KAAFLKATDE LIAAVTAHWR EDFTVLR LHG DCHAGNILWR DGPMFVDLDD ARNGPAVQDL WMLLNGDKAE QRMQLETIIE AYEEFSEFDT AEIGLIEPLR AMRLVYYLAW LMRRWADPAF PKNFPWLTGE DYWLRQTATF IEQAKVLQEP PLQLTPMY
Source	Baculovirus
Target Names	srkA
Protein Names	Recommended name: Protein rdoA
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.