



Recombinant Human Cholecystokinin (CCK), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP004771HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P06307 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | VSQRTDGESRAHLGALLARYIQQARKAPSGRMSIVKNLQNLDPShRISDRDYM GWMDF |
| Source | E.coli |
| Target Names | CCK |
| Protein Names | Recommended name: Cholecystokinin Short name= CCKCleaved into the following 10 chains: 1. Cholecystokinin-58 Short name= 2. CCK58 3. Cholecystokinin-58 desnonopeptideAlternative name(s): (1-49)-CCK58 Cholecystokinin-39 |
| Expression Region | 46-103aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Target Details | Cholecystokinin is a brain/gut peptide. In the gut, it induces the release of pancreatic enzymes and the contraction of the gallbladder. In the brain, its physiologic role is unclear. The cholecystokinin pro-hormone is processed by endo- and exo-proteolytic cleavages. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |