



Recombinant Human T-cell leukemia virus 1 Accessory protein p30II

Product Code	CSB-EP313863HQJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C214
Product Type	Recombinant Protein
Immunogen Species	Human T-cell leukemia virus 1 (strain Japan ATK-1 subtype A) (HTLV-1)
Purity	>85% (SDS-PAGE)
Sequence	MALCCFAFSA PCLHLRSRRS CSSCFLATS AAFFSARLLR RAFSSSFLFK YSAVCFSSSF SRSFFRFLFS SARRCRSRCV SPRGGAFSPG GPRRSRPRLS SSKDSKPSST ASSSSLSFNS SSKDNSPSTN SSTRSSGHD TGKHRNSPAD TKLTMLIISP LPRVWTESSF RIPSLRVWRL CTRRLVPHLW GTMFGPPTSS RPTGHLSRAS DHLGPHRWTR YRLSSTVPYP STPLLPHPEN L
Source	E.coli
Protein Names	Recommended name: Accessory protein p30II
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.