



# Recombinant Bovine DNA nucleotidylexotransferase (DNNT), partial

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP007088BO   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P06526   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Bos taurus (Bovine)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | DNNT   |
| <b>Protein Names</b>     | Recommended name: DNA nucleotidylexotransferase EC= 2.7.7.31 Alternative name(s): Terminal addition enzyme Terminal deoxynucleotidyltransferase Short name= TDT Short name= Terminal transferase   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Partial  |
| <b>Target Details</b>    | This gene is a member of the DNA polymerase type-X family and encodes a template-independent DNA polymerase that catalyzes the addition of deoxynucleotides to the 3'-hydroxyl terminus of oligonucleotide primers. In vivo, the encoded protein is expressed in a restricted population of normal and malignant pre-B and pre-T lymphocytes during early differentiation, where it generates antigen receptor diversity by synthesizing non-germ line elements (N-regions) at the junctions of rearranged Ig heavy chain and T cell receptor gene segments. Alternatively spliced transcript variants encoding different isoforms of this gene have been described. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |