



Recombinant *Oryza sativa* subsp. japonica MADS-box transcription factor 56 (MADS56)

Product Code	CSB-EP313962OFG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C5B2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MVRGRTELKR IENPTSRQVT FSKRRNGLLK KAFELSVLCD AEVALIVFSP RGRLYEFASA PSLQKTIDRY KAYTKDHVNN KTIQQDIQQV KDDTLGLAKK LEALDESRRK ILGENLEGFS IEELRGLEMK LEKSLHKIRL KKTELLEQQI AKLKEKERTL LKDNENLRGK HRNLEAAALV ANHMTTTTAP AAWPRDVPMT SSTAGAADAM DVETDLYIGL PGTERSSNRS ETG
Source	E.coli
Target Names	MADS56
Protein Names	Recommended name: MADS-box transcription factor 56 Alternative name(s): FDRMADS8 OsMADS56 RMADS214
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.