



Recombinant Human GATA-type zinc finger protein 1 (ZGLP1)

Product Code	CSB-YP314032HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C6A0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTEPQVGCVA CPRVHKEPAQ VGTPWPAKPR SHPRKRDPTA LLPRSLWPAC QESVTALCFL QETVERLGQS PAQDTPVLGP CWDPMALGTQ GRLLDRDSK DTQTRISQKG RRLQPPGTPS APPQRRPRKQ LNPCRGRTERV DPGFEGVTLK FQIKPDSSLQ IIPTYSLPCS SRSQESPADA VGGPAAHPGG TEAHSAGSEA LEPRRCASCR TQRTPLWRDA EDGTPLCNAC GIRYKKYGTR CSSCWLVP RK NVQPKRLCGR CGVSLDPIQE G
Source	Yeast
Target Names	ZGLP1
Protein Names	Recommended name: GATA-type zinc finger protein 1 Alternative name(s): GATA-like protein 1 Short name= GLP-1
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.