



Recombinant Hepatitis delta virus genotype III Large delta antigen

Product Code	CSB-EP314063HFY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C6M3
Product Type	Recombinant Protein
Immunogen Species	Hepatitis delta virus genotype III (isolate Peru-1) (HDV)
Purity	≥85% (SDS-PAGE)
Sequence	MSQTVARLTS KEREILEQW VEERKNRRKL EKDLRRANKK IKKLEDENPW LGNVVGLLRR KKDEDGAPPA KRPRQETMEV DSGPGRKPKA RGFTDQERRD HRRRKALENK KKQLAGGGKH LSQEEEEELR RLARDDDERE RRTAGPRPGG VNPMDGPPRG APGGGFVPSL QGVPESPFSR TEGEGIDIRGT QQFPWYGFTP PPPGYYWVPG C
Source	E.coli
Protein Names	Recommended name: Large delta antigen Short name= L-HDAg Alternative name(s): p27
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.