





Recombinant Dog Inactive pancreatic lipase-related protein 1 (PNLIPRP1)

Product Code	CSB-YP018263DO
Relevance	May function as inhibitor of dietary triglyceride digestion. Lacks detectable lipase activity towards triglycerides, diglycerides, phosphatidylcholine, galactolipids or cholesterol esters (in vitro).
Abbreviation	Recombinant Dog PNLIPRP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06857
Product Type	Recombinant Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	KEVCYEQIGCFSDAEPWAGTAIRPLKVLPWSPERIGTRFLLYTNKNPNNFQTLL PSDPSTIEASNFQTDKKTRFIHGFIDKGEENWLLDMCKNMFKVEEVNCICVDW KKGSTSYTQAANNVRVVGAAQVAQMLSMLSANYSYSPSQVQLIGHSLGAHVA GEAGSRTPGLGRITGLDPVEASFQGTPEEVRLDPTDADFVDVIHTDAAPLIPFL GFGTSQQMGHLDFFPNGGEEMPGCKKNALSQIVDLGDIWEGTRDFVACNHL RSYKYYSESILNPDGFASYPCASYRAFESNKCFCPCPDQGCPQMGHYADKFAV KTSDETQKYFLNTGDSSNFARWRYGVSTITLSGKRATGQAKVALFGSKGNTHQ FNIFKGILKPGSTHSNEFDAKLDVGTIEKVKFLWNNNVNPTFPKVGAAKITVQ KGEEKTVHSFCSESTVREDVLLTLTPC
Research Area	Signal Transduction
Source	Yeast
Target Names	PNLIPRP1
Protein Names	Short name:PL-RP1
Expression Region	18-467aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	51.7 kDa
Protein Length	Full Length of Mature Protein
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.