



Recombinant Mouse Transcriptional activator Myb (Myb)

Product Code	CSB-YP015261MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P06876
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MARRPRHSIY SSDEDDIEDIE MCDHDYDGLL PKSGKRHLGK TRWTREEDEK LKKLVEQNGT DDWKVIANYL PNRTDVQCQH RWQKVLNPEL IKGPWTKEED QRVIELVQKY GPKRWSVIAK HLGKGRIGKQC RERWHNHLNP EVKKTSTWTEE EDRIIYQAHK RLGNRWAEIA KLLPGRTDNA IKNHWNSTMR RKVEQEGYLQ EPSKASQTPV ATSFQKNNHL MFGGHASPPS QLSPSGQSSV NSEYPYHIA EAQNISSHVP YPVALHVNIIV NVPQAAAAI QRHYNDEDPE KEKRIKELEL LLMSTENELK GQQALPTQNH TCSYPGWHTS SIVDQTRPHG DSAPVSCSCLGE HHATPSLPAD PGSLEESAS PARCMIVHQQ TILDNVKNLL EFAETLQFID SFLNTSSNHE SSGLDAPTLP STPLIGHKLT PCRDQTVKTQ KENSIFRTPA IKRSILESSP RTPTPFKHAL AAQEIKYGPL KMLPQTPSHA VEDLQDVIKQ ESDESGIVAE FQESGPPLLK KIKQEVESPT EKSGNFFCSN HWAENSLSTQ LFSQASPVAD APNILTSSVL MTPVSEDEDN VLKAFVTPKN RPLVGPLQPC SGAWEPASCG KTEDQMTASG PARKYVNAFS ARTLVM
Source	Yeast
Target Names	Myb
Protein Names	Recommended name: Transcriptional activator Myb Alternative name(s): Proto-oncogene c-Myb
Expression Region	1-636
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a transcription factor that is a member of the MYB family of transcription factor genes. The protein contains three domains, an N-terminal DNA-binding domain, a central transcriptional activation domain and a C-terminal domain involved in transcriptional repression. This protein plays an essential role in the regulation of hematopoiesis and may play a role in tumorigenesis. Alternative splicing results in multiple transcript variants.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.