



Recombinant Catalase-peroxidase (katG)

Product Code	CSB-EP314168ENL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C8F0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	ADKKETQ NFYYPETLDL TPLRLHSPES NPWGADFDYA TRFQQLDMEA LKKDIKDLLT TSQDWWPADY GHYGPFFIRM AWHGAGTYRT YDGRGGASGG QQRFEPLNSW PDNVNLDKAR RLLWPVKKKY GSSISWGDLM VLTGNVALES MGFKTLGFAG GREDDWESDL VYWGPDNKPL ADNRDKNGKL QKPLAATQMG LIYVNPEGPG GKPDPLASAK DIREAFSRMA MDDEETVALI AGGHTFGKAH GAASPEKCIG AGPDGAPVEE QQLGWKNKCG TGNGKYTITS GLEGAWSTSP TQFTMQYLKN LYKYEWELHK SPAGAYQWKP KKAANIVQDA HDPSVLHPLM MFTTDIALKV DPEYKKITTR FLNDPKAFEQ AFARAWFKLT HRDMGPAARY LGNEVPAESF IWQDPLPAAD YTMIDGKDIK SLKEQVMDLG IPASELIKTA WASASTFRVT DYRGGNNGAR IRLQPEINWE VNEPEKLKKV LASLTSLQRE FNKKQSDGKK VSLADLIVLS GNAAIEDAAR KAGVELEIPF TPGRTDASQE QTDVASFSVL EPTADGFRNY YSKSRSHISP VESLIDKASQ LDLTVPEMTA LLGGLRVMDI NTNNSLGVF TDTPGVLDNK FFVNLLDMST RWSKADKEDT YNGFDRKTGA LKWKASSVDL IFSSNPELRA VAEVYASDDA RNKFIHDFVK SWNKVMNSDR FDLNNK
Source	E.coli
Target Names	katG
Protein Names	Recommended name: Catalase-peroxidase Short name= CP EC= 1.11.1.21 Alternative name(s): Peroxidase/catalase
Expression Region	24-736
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.