



Recombinant Human T-cell leukemia virus 1 HTLV-1 basic zipper factor (HBZ)

Product Code	CSB-MP314490HQQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C745
Product Type	Recombinant Protein
Immunogen Species	Human T-cell leukemia virus 1 (isolate Caribbea HS-35 subtype A) (HTLV-1)
Purity	>85% (SDS-PAGE)
Sequence	MVNFVSVGLF RCLPVPCPED LLVEELVDGL LSLEEELKDK EEEETVLDGL LSLEESRGR LRRGPPGGKA PPRGETHRDR QRRAEKRRK KKEREKEEEE QIAEYLKRKE EEKARRRKRA EEKAADFARR KQEEQERRER KWRQGAEKAK QHSARKEKMQ ELGVDGYTRQ LEGEVESLEA ERRRLLEKE DLMGEVNYWQ GRLEAMWLQ
Source	Mammalian cell
Target Names	HBZ
Protein Names	Recommended name: HTLV-1 basic zipper factor Short name= HBZ
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.