




Recombinant Platelet glycoprotein Ib alpha chain (GP1BA) for Homo sapiens (Human)

Product Code	CSB-EP009685HU(A4)
Abbreviation	Recombinant Human GP1BA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	BC027955
Alias	Antigen CD42b-alpha CD_antigen: CD42b
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	HPICEVSKVASHLEVNC DKRNLTALPPDLPKD TTILHLS ENLLYTFSLATLMPYT RLTQLNLDRCELTKLQVDGTL PVLGTLDSLHNQLQSLPLL GQTLPALTVLDVSF NRLTSLPLGALRGLGELQELYLKGNELKTLPPGLLTPTPKLEKLSLANNNLTELP AGLLNGLENLDTLLLQENSLYTIPKGFFGSHLLPFAFLHGNPWLCNCEILYFRR WLQDNAENVYVWKQGV DVKAMTSNVASVQCDNSDKFPVYKYPGKGCPTLG DEGDTDLYDYYPEEDTEGDKVRATRTVVKFPTKAHTTPWGLFYSWSTASLDS QMPSSLHPTQESTKEQTTFPPRWTPNFTLHMESITFSKTPKSTTEPTSPSTTS EPVPEPAPNMTTLEPTSPPTTPEPTSEAPSPPTTPEPTPIPTIATSPTILVSATSL ITPKSTFLT TTKPVSLLESTKKTIP EL DQPPKLRGVLQGHLESSRNDPFLHPDFC CLLPLGFYVLGLFWLLFASVVLILLLSWVGH
Research Area	Cardiovascular
Intended Use	preparation antibody, WB, ELISA And other possible applications
Source	E.coli
Target Names	GP1BA
Protein Names	Platelet glycoprotein Ib alpha chain
Expression Region	17-531aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	83.9 kDa
Protein Length	Partial of BC027955
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.