



Recombinant Protein tldD (tldD)

Product Code	CSB-MP314595EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGG9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSLNLVSEQL LAANGLKHQD LFAILGQLAE RRLDYGDLYF QSSYHESWVL EDRIIKDGSY NIDQGVGVRA ISGEKTGFAY ADQISLLALE QSAQAARTIV RDSGDGKVQT LGAVEHSPLY TSDVPLQSMS REEKLDILRR VDKVAREADK RVQEVTASLS GYELILVAA TDGTLAADVR PLVRLSVSVL VEEDGKRERG ASGGGGRFGY EFFLADLDGE VRADAWAKEA VRMALVNLSA VAAPAGTMPV VLGAGWPGVL LHEAVGHGLE GDFNRRGTSV FSGQVGELVA SELCTVDDG TMVDRRGSA IDDEGTPGQY NVLIENGILK GYMQDKLNAR LMGMTPTGNG RRESYAHLPMPRMTNTYMLP GKSTPQEIIIE SVEYGIYAPN FGGGQVDITS GKVFVSTSEA YLIENGKGTK PVKGATLIGS GIETMQQISM VGNDLKLDNG VGVCGKEGQS LPVGVGQPTL KVDNLTVGGT A
Source	Mammalian cell
Target Names	tldD
Protein Names	Recommended name: Protein tldD
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.