



# Recombinant Protein tldD (tldD)

<b>Product Code</b>	CSB-YP314595EGX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0AGG9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLNLVSEQL LAANGLKHQD LFAILGQLAE RRLDYGDLYF QSSYHESWVL EDRIIKDGSY NIDQGVGVRA ISGEKTGFAY ADQISLLALE QSAQAARTIV RDSGDGKVQT LGAVEHSPLY TSDPLQSMS REEKLDILRR VDKVAREADK RVQEVTASLS GYELILVAA TDGTLAADVR PLVRLSVSVL VEEDGKRERG ASGGGGRFGY EFFLADLDGE VRADAWAKEA VRMALVNLSA VAAPAGTMPV VLGAGWPGVL LHEAVGHGLE GDFNRRGTSV FSGQVGELVA SELCTVDDG TMVDRRGSVA IDDEGTPGQY NVLIENGILK GYMQDKLNAR LMGMTPTGNG RRESYAHLPMPRMTNTYMLP GKSTPQEIIIE SVEYGIYAPN FGGGQVDITS GKFVFSTSEA YLIENGKGTK PVKGATLIGS GIETMQQISM VGNDLKLDNG VGVCGKEGQS LPVGVGQPTL KVDNLTVGGT A
<b>Source</b>	Yeast
<b>Target Names</b>	tldD
<b>Protein Names</b>	Recommended name: Protein tldD
<b>Expression Region</b>	1-481
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.