



# Recombinant Trk system potassium uptake protein trkA (trkA)

<b>Product Code</b>	CSB-EP314596EGX-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0AGI9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKIIILGAGQ VGGTLAENLV GENNDITVVD TNGERLRTLQ DKFDLRVVQG HGSHPRVLRE AGADDADMLV AVTSSDETNM VACQVAYSLF NTPNRIARIR SPDYVRDADK LFHSDAVPID HLIAPEQLVI DNIYRLIEYP GALQVVNFAE GKVS LAVVKA YYGGPLIGNA LSTMHREMPH IDTRVAAIFR HDRPIRPQGS TIVEAGDEVF FIAASQHIRA VMSELQRLEK PYKRIMLVGG GNIGAGLARR LEKDYSVKLI ERNQQRAAEL AEKLQNTIVF FGDASDQELL AEEHIDQVDL FIAVTNDDEA NIMSAMLAKR MGAKKVMVLI QRRAYVDLVQ GSVIDIAISP QQATISALLS HVRKADIVGV SSLRRGVAEA IEAVAHGDES TSRVVGRVID EIKLPPGTII GAVVRGNDVM IANDNLRIEQ GDHVIMFLTD KKFITDVERL FQPSPFFL
<b>Source</b>	E.coli
<b>Target Names</b>	trkA
<b>Protein Names</b>	Recommended name: Trk system potassium uptake protein trkA Short name=K(+)-uptake protein trkA
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.