



Recombinant Arabidopsis thaliana Nudix hydrolase 24, chloroplastic (NUDT24)

Product Code	CSB-MP314604DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C026
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SSSRSPPLHI SRVLSVETVP LSPSFTWNDV FENSRKEYVP QNSSDLTGFL EKVDRCNRGL EKLGEFIPFV IEEQIVGYIH KGFTKYLRDF NDIFTSQYG GHVTLNMLLD KPEERTRA VA HVIKILGNKG IIPGIRNELY PVKPSFNAPA FFSIERAAAP YFGLKGYAIH VNGYVERDGQ KFLWIGKRSL AKSTYPGKLD HLVAGGLPHG ISVCENLVKE CEEEAGISKV LADRAIAGV VSYMDIDRYC FTRDVLFCYD LELPQDFVPT NQDGEVDSFR LIPVAQVANV VRKTSFFKDS CSLVIIDFLF RHGLIRPESP GYLDLYRRLR NGDCS
Source	Mammalian cell
Target Names	NUDT24
Protein Names	Recommended name: Nudix hydrolase 24, chloroplastic Short name= AtNUDT24 EC= 3.6.1.-
Expression Region	31-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.