



Recombinant Clostridium perfringens Phospholipase C (plc)

Product Code	CSB-BP314672CMB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C216
Product Type	Recombinant Protein
Immunogen Species	Clostridium perfringens
Purity	≥85% (SDS-PAGE)
Sequence	WD GKIDGTGTHA MIVTQGVSI L ENDSLKNEPE SVRKNLEILK ENMHELQLGS TYPDYDKNAY DLYQDHFWD P DTDNNFSKDN SWYLAYSIPD TGESQIRKFS ALARYEWQRG NYKQATFYL G EAMHYFGDID TPYHPANVTA VDSAGHVKFE TFAEERKEQY KINTAGCKT N EDFYADILKN KDFNAWSKEY ARGFAKTGKS IYYSHASMSH SWDDWDYAAK VTLANSQKGT AGYIYRFLHD VSEGNDPSVG KNVKELVAYI STSGEKDAGT DDYMYFGIKT KDGKTQEWEM DNPGNDFMTG SKDTYTFK LK DENLKIDDIQ NMWIRKRKYT AFPDAYKPEN IKIIANGKVV VDKDINEWIS GNSTYNIK
Source	Baculovirus
Target Names	plc
Protein Names	Recommended name: Phospholipase C Short name= PLC EC= 3.1.4.3Alternative name(s): Alpha-toxin Hemolysin Lecithinase Phosphatidylcholine cholinephosphohydrolase
Expression Region	29-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.