



Recombinant *Oryza sativa* subsp. *indica* Uncharacterized protein ycf68 (ycf68-1)

Product Code	CSB-EP314707OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C305
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAYSSCLNRS LKPNKLLRR IDGAIQVRSH VDLTFYSLVG SGRSGGGTTA PLFSRIHTSL ISVWRAISRA QVEVRPQWEN GAPNNASSQT KNYEITLSFW GDGGIVPFEP FFHAFPGGLE KAAINRTSLI LPS
Source	E.coli
Target Names	ycf68-1;
Protein Names	Recommended name: Uncharacterized protein ycf68
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.