



Recombinant Ribosome association toxin RatA (ratA)

Product Code	CSB-MP314729SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGL6
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MILFVGFLLM EIVMPQISRT ALVPYSAEQM YQLVNDVQSY PQFLPGCTGS RILESTPGQM TAAVDVSKAG ISKTFTTRNQ LTSNQSILMN LVDGPFKKLI GGWKFTPLSQ EACRIEFHLD FEFTNKLIEL AFGRVFKELA ANMVQAFTVR AKEVYSAR
Source	Mammalian cell
Target Names	ratA
Protein Names	Recommended name: Ribosome association toxin RatA
Expression Region	1-158
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.