



# Recombinant Clostridium perfringens Perfringolysin O (pfo)

<b>Product Code</b>	CSB-MP314820CMB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0C2E9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium perfringens (strain 13 / Type A)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	KD ITDKNQSIDG GISSLSYNRN EVLASNGDKI ESFVPKEGKK TGNKFIVVER QKRSLTTSPV DISIIDSVND RTYPGALQLA DKAFVENRPT ILMVKRKPIN INIDLPLGKG ENSIKVDDPT YGKVSGAIDE LVSKWNEKYS SHTLTPARTQ YSESMVYSKS QISSALNVNA KVLNSLQVD FNAVANNEKK VMILAYKQIF YTVSADLPKN PSDLFDDSVT FNDLKQKQGVN NEAPPLMVSN VAYGRTIYVK LETTSSSKDV QAAFKALIKN TDIKNSQQYK DIYENSSFTA VVLGGDAQEH NKVVTKDFDE IRKVIKDNAT FSTKNPAYPI SYTSVFLKDN SVAAVHNKTD YIETTSTEYS KGKINLDHSG AYVAQFEVAW DEVSYDKEGN EVLTHKTWDG NYQDKTAHYS TVIPLEANAR NIRIKARECT GLAWEWWRDV ISEYDVPLTN NINVSIWGTT LYPGSSITYN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pfo
<b>Protein Names</b>	Recommended name: Perfringolysin O Alternative name(s): Theta-toxin Thiol-activated cytolysin
<b>Expression Region</b>	29-500
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.