



Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At1g06630 (At1g06630)

Product Code	CSB-MP314822DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C2G0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDAKKLDTGP RDAINWLPDE ILGKILSLLA TKQAVSTSVL SKKWRTLFLKLVDTLEFDDSV SGMGEQEASY VFPESFKDLV DRTVALQCDY PIRKLSLKCHVGRDDEQRKA CVGRWISNVV GRGVSEVVL R INDRGLHFLS PQLLTCKTLV KLTLGTRLFL GKLPYSVSLP SLKFLFIHSV FFDDFGELSN VLLAGCPVVE ALYLNQNGES MPYTISSPTL KRLSVHYEYH FESVISFDLP NLEYLDYSDY ALYGYPQVNL ESLVEAYLNL DKAEHVESPD VTKLIMGIRN VEILSLSPDS VGVYISCKY GLLLPVFNNL VSLSFGTKKT RAWKLLADIL KQSPKLETLI IEDLNGYPLD VSMPLNQVKE LQILEYGESD DEVKRLKSFL GEESMIEVVF PEV
Source	Mammalian cell
Target Names	At1g06630
Protein Names	Recommended name: F-box/LRR-repeat protein At1g06630
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.