



Recombinant *Oryza sativa* subsp. japonica Auxin-responsive protein IAA31 (IAA31)

Product Code	CSB-BP314909OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C133
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MENLKATELR LGLPGTEEEA APPPSTPRAG SKRALAGEPD QAKIKPAAAA KAQVVGWPPV RSYRKSLQP TTTTTSKPP PAAAAAETQQ KEDVAGAGGL FVKVSMGDGAP YLRKIDLVY KGYRELREAL EAMFLCFSGG AAADAAVNPS DFAVTYEDKD GDLMLVGDVP FEMFISTCKR LRIMKGSEAR GLGATRG
Source	Baculovirus
Target Names	IAA31
Protein Names	Recommended name: Auxin-responsive protein IAA31 Alternative name(s): Indoleacetic acid-induced protein 31
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.