



Recombinant Zea mays Expansin-B11 (EXPB11)

Product Code	CSB-MP314944ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C1Y5
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	GPPKVPPGKN ITAKYGSDWL DAKATWYGKP TGAGPDDNGG GCGYKDVNKA PFNSMGACGN VPIFKDGLGC GSCFEIKCDK PAECSGKPVV VYITDMNYEP IAAYHFDLAG TAFGAMAKKG EEEKLRKAGI IDMQFRRVKC KYGSKVTFHL EKGCNPNYLA LLVKYVDGDG DIVAVDIKEK GSDTYEPLKH SWGAIWRKDS DKPIKGPITV QLTTEGGTKT VYDDVIPAGW KPNTAYTAK
Source	Mammalian cell
Target Names	EXPB11
Protein Names	Recommended name: Expansin-B11 Alternative name(s): Beta-expansin-11 Pollen allergen Zea m 1a Pollen allergen Zea m 1b ZmEXPB11 Allergen= Zea m 1
Expression Region	31-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.