



# Recombinant Rat 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 1 (Pfkfb1)

<b>Product Code</b>	CSB-MP017817RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P07953
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SREMGELTQ TRLQKIWIPH SSSSSVLQRR RGSSIPQFTN SPTMVIMVGL PARGKTYIST KLTRYLNWIG TPTKVFNLGQ YRREAVSYRN YEFFRPDNT AQLIRKQCAL AALKDVHKYL SREEGHVAVF DATNTTRERR SLILQFAKEH GYKVFFIESI CNDPEIIAEN IKQVKLGSPD YIDCDQEKVL EDFLKRIECY EINYQPLDEE LDSHLSYIKI FDVGTRYMVN RVQDHSVQSR AYYLMNIHVT PRSIYLCRHG ESELNLRGRI GGDSGLSARG KQYAYALANF IRSQGISSLK VWTSHMKRTI QTAEALGVPY EQWKALNEID AGVCEEMTYE EIQEHYPEEF ALRDQDKYRY RYPKGESYED LVQRLEPVIM ELERQENVLV ICHQAVMRC LAYFLDKSSD ELPYLKCPH TVLKLTPVAY GCRVESIYLN VEAVENTHRDK PENVDITREA EEALDTVPAH Y
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Pfkfb1
<b>Protein Names</b>	Recommended name: 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 1 Short name= 6PF-2-K/Fru-2,6-P2ase 1 Short name= PFK/FBPase 1 Alternative name(s): 6PF-2-K/Fru-2,6-P2ase liver isozyme Including the following 2 domains:
<b>Expression Region</b>	2-471
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of the family of bifunctional 6-phosphofructo-2- kinase:fructose-2,6-bisphosphatase enzymes. The enzyme forms a homodimer that catalyzes both the synthesis and degradation of fructose-2,6-bisphosphate using independent catalytic domains. Fructose-2,6-bisphosphate is an activator of the glycolysis pathway and an inhibitor of the gluconeogenesis pathway. Consequently, regulating fructose-2,6-bisphosphate levels through the activity of this enzyme is thought to regulate glucose homeostasis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.